



Mouse Anti-Chicken Collagen type II Monoclonal antibody, DJJD2 (CABT-S-H2079)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	native chick type I1 collagen (Full length protein)
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Bovine, Chicken, Fish, Human, Mouse, Ovine, Rabbit, Rat, Xenopus
Clone	DJJD2
Purification	Unpurified
Conjugate	Unconjugated
Applications	ELISA, IF, IHC, WB
Format	Liquid
Concentration	Lot specific
Size	1 ml
Preservative	None
Storage	Maintain at -20°C for up to 12 months. Avoid repeated freeze-thaw cycles. Store product undiluted.
Ship	Wet ice

BACKGROUND

Introduction	Type II collagen is specific for cartilaginous tissues. It is essential for the normal embryonic development of the skeleton, for linear growth and for the ability of cartilage to resist compressive forces
Keywords	COL2A1; collagen, type II, alpha 1; collagen alpha-1(II) chain; chondrocalcin; alpha-1 type II collagen; disproportionate micromelia