



Rat Anti-Human IL6 monoclonal antibody, clone NR3-24B6 (CABT-WN1691)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Monoclonal Antibody to Human IL6 molecule
Antigen Description	IL-6 is a potent lymphoid cell growth factor that stimulates the growth and survival of certain B cells and T cells. IL-6 plays a role in host defense, acute phase reactions, immune response, and hematopoiesis. IL-6 is expressed by T cells, B cells, monocytes, fibroblasts, hepatocytes, endothelial cells, and keratinocytes.
Target	IL6
Immunogen	COS-7- expressed, recombinant human IL-6
Isotype	IgG1
Source/Host	Rat
Species Reactivity	Human
Clone	NR3-24B6
Purification	The LEAF (Low Endotoxin, Azide-Free) antibody was purified by affinity chromatography.
Conjugate	Functional Grade
Applications	Neut, ICFC, ELISA(Cap), ELISPOT
Format	Liquid
Concentration	0.5 mg/ml
Size	500 µg

Buffer	PBS only.
--------	-----------

Preservative	None
--------------	------

GENE INFORMATION

Entrez Gene ID	3569
----------------	----------------------

UniProt ID	B4DVM1
------------	------------------------
