



# Rat anti-Human IL6 Monoclonal antibody, clone NR3-24B6 (CABT-WN1697)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rat Monoclonal Antibody to Human IL6 molecule
<b>Antigen Description</b>	IL-6 is a potent lymphoid cell growth factor that stimulates the growth and survival of certain B cells and T cells. IL-6 plays a role in host defense, acute phase reactions, immune response, and hematopoiesis. IL-6 is expressed by T cells, B cells, monocytes, fibroblasts, hepatocytes, endothelial cells, and keratinocytes.
<b>Target</b>	IL6
<b>Immunogen</b>	COS-7- expressed, recombinant human IL-6
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Rat
<b>Species Reactivity</b>	Human
<b>Clone</b>	NR3-24B6
<b>Purification</b>	The antibody was purified by affinity chromatography.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	CyTOF , ELISA(Cap), ELISPOT
<b>Format</b>	Liquid
<b>Concentration</b>	0.2 mg/ml
<b>Size</b>	500 µg

<b>Buffer</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
---------------	---

<b>Preservative</b>	0.09% Sodium Azide
---------------------	--------------------

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">3569</a>
-----------------------	----------------------

<b>UniProt ID</b>	<a href="#">B4DVM1</a>
-------------------	------------------------