



Rabbit Anti-Guinea Pig Urinary Protein Polyclonal Antibody (CABT-YN1491)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Binds to guinea pig urinary proteins. Reacts with guinea pig epithelial and hair extracts. Does not recognize rat or mouse urinary proteins; does not react with rabbit or hamster epithelial and hair extracts.
Target	Guinea Pig Urinary Protein
Immunogen	Guinea pig urinary protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Guinea pig
Purification	Affinity chromatography using recombinant Protein G.
Conjugate	Unconjugated
Applications	ELISA Recommended dilution: ELISA: 1:1000 Final working dilutions must be determined by end user.
Format	Liquid
Concentration	Lot specific
Size	100 µl
Preservative	None

Storage	Store at 4°C
----------------	--------------

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Allergic reaction to guinea pig has been recognized as a problem in domestic settings and work environments for many years. Three major allergens were identified within both sources with molecular weights (MW) of 8 kDa, 17 kDa and 20 kDa, respectively.
---------------------	--

Keywords	guinea pig urinary protein; Urinary Protein
-----------------	---

GENE INFORMATION

Protein Refseq	None
-----------------------	------