



Mouse Anti-Perennial Ryegrass Lol p 1 Monoclonal Antibody, clone 8C10 (CABT- YN1521)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Binds to a common epitope on rye grass pollen <i>Lolium perenne</i> allergen Lol p 1
Target	Perennial Ryegrass Lol p 1
Immunogen	Lol p 1
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	<i>Lolium perenne</i>
Clone	8C10
Purification	Produced in tissue culture and purified by affinity chromatography using Protein G.
Conjugate	Unconjugated
Applications	ELISA Recommended dilution: ELISA: 1:1000 Final working dilutions must be determined by end user.
Format	Liquid
Concentration	Lot specific
Size	100 μ l

Buffer	Phosphate buffered saline, pH 7.4, preservative free.
Preservative	None
Storage	Store at 4°C
Ship	Wet ice

BACKGROUND

Introduction	Lol p 1 is a group 1 pollen allergen (beta expansin).
Keywords	Lol p1; Lol p 1; Lolium perenne

GENE INFORMATION

Protein Refseq	None
-----------------------	------