



Rabbit anti-Arabidopsis thaliana EIN3 (C-term) Polyclonal Antibody (CABT-Z063R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of arabidopsis thaliana EIN3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:500-1:2,000 (detect endogenous protein*)
Molecular Weight	Predicted M.W.: 71 kDa; Observed M.W.: 90 kDa
Preparation	Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of arabidopsis thaliana EIN3 (At3g20770).
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	Supplied in 1 x PBS (pH 7.4), 100 ug/ml BSA, 40% Glycerol, 0.01% NaN ₃ .
Preservative	0.01% NaN ₃

Storage	Store at -20°C. Stable for 6 months from date of receipt.
----------------	---

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	EIN3 is a probable transcription factor acting as a positive regulator in the ethylene response pathway. EIN3 is required for ethylene responsiveness in adult plant tissues. EIN3 binds to a primary ethylene response element present in the ETHYLENE-RESPONSE-FACTOR1 promoter with consequence to activate the transcription of this gene.
---------------------	--

Keywords	Protein ETHYLENE INSENSITIVE 3;At3g20770;MOE17.8;ATEIN3
-----------------	---

GENE INFORMATION

Gene Name	EIN3
------------------	------

Entrez Gene ID	821625
-----------------------	------------------------

UniProt ID	O24606
-------------------	------------------------