



Rabbit Anti-Arabidopsis thaliana eEF1B-alpha 2 Polyclonal Antibody (CABT-Z1005R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH-conjugated synthetic peptide of eEF1B-alpha 2 derived from Arabidopsis thaliana AT5G19510.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Unpurified
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:1000-1:2000
Molecular Weight	24 kDa
Reconstitution	Reconstitution with 150 µl of sterile water.
Format	Lyophilized
Size	150 µg
Preservative	None
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.

Ship

Wet ice

BACKGROUND

Introduction

eEF1B-alpha 2 is one of Translation elongation factor EF1B/ribosomal protein S6 family protein. It is the primary nucleotide exchange factor.

Keywords

eEF1B-alpha 2;Elongation factor 1-beta 2;EF-1-beta 2

GENE INFORMATION

Gene Name

At5g19510

Entrez Gene ID

[832071](#)

UniProt ID

[Q9SCX3](#)
