



Rabbit Anti-Arabidopsis thaliana FD3 Polyclonal Antibody (CABT-Z1026R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH-conjugated synthetic peptide of FD3 derived from Arabidopsis thaliana AT2G27510.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Unpurified
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:1000-1:2000
Molecular Weight	Expected M.W.: 17 kDa; Apparent M.W.: 11 kDa
Reconstitution	Reconstitution with 150 µl of sterile water.
Cellular Localization	Chloroplast
Format	Lyophilized
Size	150 µg
Preservative	None
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.

Ship

Wet ice

BACKGROUND

Introduction

Ferredoxin 3 (FD3) localized in chloroplast, it is involved in electron transport chain.

Keywords

FD3;ATFD3;FERREDOXIN 3

GENE INFORMATION

Gene Name

FD3

Entrez Gene ID

[817297](#)

UniProt ID

[Q9ZQG8](#)
