



# Rabbit Anti-Arabidopsis thaliana HEMA2 Polyclonal Antibody (CABT-Z1092R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	KLH-conjugated synthetic peptide of HEMA2 derived from Arabidopsis thaliana AT1G09940.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Purification</b>	Unpurified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB Recommended dilution: WB: 1:1000-1:2000
<b>Molecular Weight</b>	Expected M.W.: 58 kDa; Apparent M.W.: 52 kDa
<b>Reconstitution</b>	Reconstitution with 150 µl of sterile water.
<b>Cellular Localization</b>	Chloroplast
<b>Format</b>	Lyophilized
<b>Size</b>	150 µg
<b>Preservative</b>	None
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.

## BACKGROUND

**Introduction**

Glutamyl-tRNA reductase is involved in heme biosynthesis in non-photosynthetic tissues and induced by oxidative stress in photosynthetic tissues to supply heme for defensive hemoproteins.

---

**Keywords**

HEMA2;GluTR

---

## GENE INFORMATION

**Gene Name**

HEMA2

---

**Entrez Gene ID**

[837528](#)

---

**UniProt ID**

[P49294](#)

---