



Rabbit anti-Arabidopsis thaliana PHYB (C-term) Polyclonal Antibody (CABT-Z109R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of arabidopsis thaliana PHYB.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:500-1:2,000 (detect endogenous protein*)
Molecular Weight	129 kDa
Preparation	Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of arabidopsis thaliana PHYB (At2g18790).
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	Supplied in 1 x PBS (pH 7.4), 100 ug/ml BSA, 40% Glycerol, 0.01% NaN ₃ .
Preservative	0.01% NaN ₃

Storage	Store at -20°C. Stable for 6 months from date of receipt.
----------------	---

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Phytochrome B (PHYB) is a regulatory photoreceptor, which exists in two forms that are reversibly interconvertible by light: the Pr form that absorbs maximally in the red region of the spectrum and the Pfr form that absorbs maximally in the far-red region. Photoconversion of Pr to Pfr induces an array of morphogenetic responses, whereas reconversion of Pfr to Pr cancels the induction of those responses. Pfr controls the expression of a number of nuclear genes including those encoding the small subunit of ribulose-bisphosphate carboxylase, chlorophyll A/B binding protein, protochlorophyllide reductase, rRNA, etc. It also controls the expression of its own gene(s) in a negative feedback fashion. PHYA is involved in the flowering time regulation.
---------------------	---

Keywords	Phytochrome B;HY3
-----------------	-------------------

GENE INFORMATION

Gene Name	PHYB
------------------	------

Entrez Gene ID	816394
-----------------------	------------------------

UniProt ID	P14713
-------------------	------------------------