



Rabbit anti-Arabidopsis thaliana TMK1 (Middle region) Polyclonal Antibody (CABT-Z145R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antibodies were produced by immunizing animals with a GST-fusion protein containing the middle region of arabidopsis thaliana TMK1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:500-1:2,000 (detect endogenous protein*)
Molecular Weight	102 kDa
Preparation	Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the middle region of arabidopsis thaliana TMK1 (At1g66150).
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	Supplied in 1 x PBS (pH 7.4), 100 ug/ml BSA, 40% Glycerol, 0.01% NaN ₃ .
Preservative	0.01% NaN ₃

Storage	Store at -20°C. Stable for 6 months from date of receipt.
----------------	---

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	TRANSMEMBRANE KINASE 1 (TMK1) is a receptor-like transmembrane protein serine/threonine kinase. TMK1 is involved in protein phosphorylation and signal transduction. TMK1 has ATP binding and protein serine/threonine kinase activity.
---------------------	---

Keywords	Probable receptor protein kinase TMK1;TRANSMEMBRANE KINASE 1;BARK1-like kinase 1;Transmembrane kinase 1
-----------------	---

GENE INFORMATION

Gene Name	TMK1
------------------	------

Entrez Gene ID	842930
-----------------------	------------------------

UniProt ID	P43298
-------------------	------------------------