



Rabbit Anti-Human CDH1 Monoclonal Antibody, clone CQ7159 (CABT-Z219R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to E-Cadherin residues within aa600-700 of E-cadherin was used as an immunogen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	CQ7159
Purification	ProA affinity purified IgG.
Conjugate	Unconjugated
Applications	IHC-P Recommended concentration: IHC-P: 1:100-1:200
Molecular Weight	97 kDa
Cellular Localization	Membrane
Positive Control	Liver
Format	Liquid
Concentration	Lot specific
Size	100 µl

Buffer	PBS 59%, Sodium azide 0.01%, Glycerol 40%, BSA 0.05%.
Preservative	0.01% Sodium azide
Storage	Store at -20 °C. Avoid freeze/thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	E-cadherin is a transmembrane, calcium dependent cell adhesion protein that mediates cell to cell adhesion and maintains structural and functional integrity of epithelial tissues. It also has pivotal barrier functions and maintains the polarity of epithelial cells. Reduced or aberrant Ecadherin expression breaks cell to cell contacts, and thus, cells acquire the ability to migrate. In normal tissues, immunostaining of E-cadherin is localized to the membrane of epithelial cells, consistent with its role in cell adhesion. And in tumor tissues, E-cadherin stains positively in glandular epithelium as well as adenocarcinomas of the lung, gastrointestinal tract, and ovary. It has also been shown to be positive in some thyroid carcinomas. A combination of E-cadherin and p120 catenin may help distinguish ductal carcinoma of the breast from lobular carcinoma. And also, several studies have reported that reduced E-cadherin expression is correlated with poor prognosis in several types of carcinomas.
Keywords	CDH1; cadherin 1, type 1, E-cadherin (epithelial); UVO; CDHE; ECAD; LCAM; Arc1; CD324; cadherin-1; CAM 120/80; E-Cadherin; uvomorulin; cell-CAM 120/80; epithelial cadherin; cadherin 1, E-cadherin (epithelial); calcium-dependent adhesion protein, epithelial

GENE INFORMATION

Gene Name	CDH1
Entrez Gene ID	999
UniProt ID	P12830