



Rabbit Anti-Human CD45 Monoclonal Antibody, clone CQ7184 (CABT-Z240R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment within Human CD45 aa-600.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	CQ7184
Purification	ProA affinity purified IgG.
Conjugate	Unconjugated
Applications	IHC-P Recommended concentration: IHC-P: 1:100-1:200
Molecular Weight	147 kDa
Cellular Localization	Membrane
Positive Control	Tonsil
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	PBS 59%, Sodium azide 0.01%, Glycerol 40%, BSA 0.05%.

Preservative	0.01% Sodium azide
Storage	Store at -20 °C. Avoid freeze/thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	CD45 is a family of single chain transmembrane glycoproteins, consisting of at least four isoforms which share a common large intracellular domain. Their extracellular domains are rod shaped, heavily glycosylated. The different isoforms are produced by alternative messenger RNA splicing of three exons of a single gene on chromosome 1q31-32. CD45 is exclusively expressed in haematolymphoid cells. Almost all haematolymphoid cells, including precursor cells and mature B- and T-lymphocytes, granulocytes, monocytes/histiocytes and interdigitating reticulum cells and follicular dendritic cells express CD45. This protein is detected in the large majority of haematolymphoid neoplasms, i.e., leukemias and malignant lymphomas. Overall, about 90% of malignant lymphomas are CD45 positive. CD45 is an important marker in the primary tumor screening panel in order to identify haematolymphoid differentiation. Loss of CD45 in precursor B-cell neoplasms is a negative prognostic parameter.
---------------------	--

Keywords	CD45; Leukocyte Common Antigen; LCA
-----------------	-------------------------------------

GENE INFORMATION

Gene Name	CD45
Entrez Gene ID	5788
UniProt ID	P08575