



Rabbit Anti-Human IL2RA Monoclonal Antibody, clone CQ7231 (CABT-Z275R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment corresponding to Human IL-2 Receptor alpha with aa1-250 was used as an immunogen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	CQ7231
Purification	ProA affinity purified IgG.
Conjugate	Unconjugated
Applications	IHC-P Recommended concentration: IHC-P: 1:100-1:200
Molecular Weight	31 kDa
Cellular Localization	Membrane
Positive Control	Colon Carcinoma Tissue
Format	Liquid
Concentration	Lot specific
Size	100 µl

Buffer	PBS 59%, Sodium azide 0.01%, Glycerol 40%, BSA 0.05%.
Preservative	0.01% Sodium azide
Storage	Store at -20 °C. Avoid freeze/thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	CD25, or interleukin 2 (IL2) receptor alpha (IL2RA), and beta (IL2RB) chains, together with the common gamma chain (IL2RG), constitute the high-affinity IL2 receptor. Homodimeric alpha chains (CD25) result in low-affinity receptor, while homodimeric beta (IL2RB) chains produce a medium-affinity receptor. Normally an integral-membrane protein, soluble CD25 has been isolated and determined to result from extracellular proteolysis. CD25 (the p55 chain of the interleukin-2 receptor) is a marker of activation for a number of cell types, including T cells, B cells, and macrophages. It also appears to be a reliable immunohistochemical marker for the discrimination of neoplastic from normal/reactive mast cells, with potential as a diagnostic tool in systemic mastocytosis. CD25 is also expressed on hairy cell leukemia.
Keywords	IL2RA; CD25; interleukin-2 receptor subunit alpha; p55; IL2-RA; IL-2-RA; IL-2R alpha; TAC antigen; cytokine receptor; IL-2R subunit alpha; IL-2 receptor alpha subunit; IL-2 receptor subunit alpha; interleukin-2 receptor alpha chain

GENE INFORMATION

Gene Name	IL2RA
Entrez Gene ID	3559
UniProt ID	P01589