



Rabbit Anti-Mouse CLDN13 Polyclonal Antibody (CABT-Z399R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide (14-aa) derived from the C-terminal intracellular region of mouse Claudin-13 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IF
Format	Liquid
Concentration	Lot specific
Size	200 µl
Buffer	PBS, pH 7.2, 0.1% sodium azide.
Preservative	0.1% sodium azide
Storage	-20°C or below.
Ship	Wet ice

BACKGROUND

Introduction

Claudin 13 is a member of claudin family. It is a four-transmembrane protein with WWCC motif. This protein plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.

Keywords

Claudin-13;Cldn13

GENE INFORMATION

Gene Name

Cldn13

Entrez Gene ID

[57255](#)

UniProt ID

[Q9Z0S4](#)