



# Rabbit Anti-Mouse CLDN13 Polyclonal Antibody (CABT-Z400R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Immunogen** Synthetic peptide directed towards the C terminal of mouse CLDN13 (NP\_065250). Peptide sequence  
SSLSLAWTSSLLLLLGGILLCVNIPVCRDFPRCIETPSARPSGANNDTLD.  
The peptide sequence for this immunogen was taken from within the described region.

**Isotype** IgG

**Source/Host** Rabbit

**Species Reactivity** Mouse, Human

**Purification** Protein A purified

**Conjugate** Unconjugated

**Applications** WB

**Format** Liquid

**Concentration** Lot specific

**Size** 100 µl

**Buffer** PBS and 2% Sucrose.

**Preservative** 0.09% sodium azide

**Storage** Store at -20°C long term. Avoid freeze / thaw cycle.

**Ship** Wet ice

# BACKGROUND

**Introduction** Claudin 13 is a member of claudin family. It is a four-transmembrane protein with WWCC motif. This protein plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.

**Keywords** Claudin-13;Cldn13

# GENE INFORMATION

**Gene Name** Cldn13

**Entrez Gene ID** [57255](#)

**UniProt ID** [Q9Z0S4](#)