



# Sheep Anti-N-Carboxymethyl LSD Polyclonal Antibody (CABT-Z411S)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	N-Carboxymethyl LSD-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

## BACKGROUND

<b>Introduction</b>	LSD is a crystalline compound, C <sub>20</sub> H <sub>25</sub> N <sub>3</sub> O, derived from lysergic acid and used as a powerful hallucinogenic drug.
<b>Keywords</b>	Carboxymethyl LSD;LSD;Lysergic acid;Lysergic acid diethylamide