



Sheep Anti-AHD Polyclonal Antibody (CABT-Z413S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Aminohydantoin-BSA
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction	1-Aminohydantoin hydrochloride may be used as one of the reactants in the synthesis of (E)-1-(2 hydroxybenzylideneamino)imidazolidine-2,4-dione, dantrolene and dantrolene sodium.
---------------------	--

Keywords

1-Aminohydantoin;AHD;1-Aminohydantoin hydrochloride;2827-56-7;1-aminoimidazolidine-2,4-dione hydrochloride;1-Aminohydantoin HCl;1-Amino hydantoin hydrochloride;2,4-Imidazolidinedione, 1-amino-, monohydrochloride;AC1Q3DFU;ACMC-209h1m;C3H6ClN3O2;1-aminohydantoin hydrochloride
