



# Sheep Anti-Chlormadinone Acetate Polyclonal Antibody (CABT-Z414S)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Chlormadinone Acetate-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

## BACKGROUND

Introduction	Chlormadinone is a steroidal progestin which was never marketed. An acylated derivative, chlormadinone acetate, is widely used clinically used as a pharmaceutical drug. While
--------------	--

chlormadinone is sometimes used as a synonym for chlormadinone acetate, what is almost always being referred to is actually chlormadinone acetate and not chlormadinone.

---

<b>Keywords</b>	Chlormadinone;Chlormadinone acetate
-----------------	-------------------------------------

---