



Sheep Anti-Codeine Polyclonal Antibody (CABT-Z417S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Dehydromethyltestosterone-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction	Dehydromethyltestosterone is a synthetic anabolic steroid with low androgenic effects and little pregestational activity. Dehydromethyltestosterone has increased in popularity amongst
--------------	---

bodybuilders.

Keywords

Nerobolettes;NSC-42722;Protobolin;Stenolon;Stenolone;Sterolon;Testosterone, 1-dehydro-17-methyl-;Metandiemone;1, 2-Dehydro-17-methyltestosterone;17-alpha-methyl-17-beta-hydroxy-1, 4-androstadien-3-one
