



Sheep Anti-Desmycosin Polyclonal Antibody (CABT-Z420S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Desmycosin-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction	Desmycosin is a 16-membered macrolide antibiotic and a metabolite of Tylosin. Desmycosin is also a precursor to Tilmicosin.
--------------	---

Keywords

Desmycosin;C39H65NO14;Tylosin B
