



# Sheep Anti-Dicloxicillin Polyclonal Antibody (CABT-Z423S)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Dicloxicillin-BSA
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

## BACKGROUND

Introduction	Dicloxacillin is a prescription drug licensed to treat infections caused by certain types of bacteria.
--------------	--

**Keywords**Dicloxicillin

---