



Sheep Anti-Fumonisin Polyclonal Antibody (CABT-Z432S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fumonisin-KLH
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction	Fumonisin are mycotoxins from <i>Fusarium moniliforme</i> that are common grain contaminants responsible for food poisoning episodes all over the world.
--------------	--

Keywords

Fumonisin
