



Sheep Anti-Methyltestosterone Polyclonal Antibody (CABT-Z445S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Methyltestosterone-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Methyltestosterone, sold under the brand names Android, Metandren, and Testred among others, is an androgen and anabolic steroid (AAS) medication which is used in the treatment of low testosterone levels in men, delayed puberty in boys, at low doses as a component of menopausal hormone therapy for menopausal symptoms like hot flashes, osteoporosis, and low sexual desire in women, and to treat breast cancer in women.

Keywords

Testoviron;methyltestosterone