



Sheep Anti-Procaïnamide Polyclonal Antibody (CABT-Z450S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Procaïnamide-BSA
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Procainamide is a pharmaceutical antiarrhythmic agent used for the medical treatment of cardiac arrhythmias, classified by the Vaughan Williams classification system as class Ia.

Keywords

Pronestyl;Procan;Procanbid;Procainamide;4-amino-N-(2-diethylaminoethyl) benzamide
