



Sheep Anti-Salicylate Polyclonal Antibody (CABT-Z453S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Salicylic Acid-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction	Salicylate is a monohydroxybenzoate that is the conjugate base of salicylic acid. It has a role as a plant metabolite. It is a conjugate base of a salicylic acid.
Keywords	salicylate;2-Hydroxybenzoate;o-Hydroxybenzoate;Salicylate ion;Salicylate anion;hydroxybenzoate