



Sheep Anti-Stilbene Monoclonal Antibody, clone 21414 (CABT-Z467S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Stilbene-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Clone	21414
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	100 µg, 1 mg
Buffer	PBS with sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Stilbenes are low-molecular weight compounds that are found in a wide range of natural sources and that exhibit a broad spectrum of biological activities, as well as application in molecular photonics and optoelectronics.

Keywords

(E)-Stilbene; 1,1'-[(E)-Ethene-1,2-diyl]dibenzene; Bibenzylidene; trans- α,β -Diphenylethylene; (E)-1,2-Diphenylethylene; ((1E)-2-Phenylvinyl)benzene; trans-Stilbene; [(E)-2-Phenylethenyl]benzene