



Mouse Anti-Enramycin Monoclonal Antibody, clone ERM (CABT-Z526M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Enramycin-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	ERM
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Enramycin or Enduracidin is a polypeptide antibiotic produced by *Streptomyces fungicidus*. Enramycin is widely used as a feed additive for pigs and chickens to prevent necrotic enteritis induced by Gram-positive gut pathogens.

Keywords

Enramycin;ERM;Enduracidin;C107H138Cl2N26O31
