



Mouse Anti-Acid Orange Monoclonal Antibody, clone AO1 (CABT-Z529M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Acid Orange-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	AO1
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Acid orange 1 is an organic compound and an azo dye.

Keywords

Acid Orange1;C.I.13090;Akacid Milling;Yellow G
