



Mouse Anti-Tebuconazole Monoclonal Antibody, clone TBZ (CABT-Z530M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Tebuconazole-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	TBZ
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Tebuconazole is a triazole fungicide used agriculturally to treat plant pathogenic fungi.

Keywords

Tebuconazole;(RS)- 1-(4-Chlorophenyl)- 4,4-dimethyl-3-(1H, 1,2,4-triazol-1-ylmethyl)pentan- 3-ol;(±)-1-(4-Chlorophenyl)-4,4-dimethyl-3-(1H, 1,2,4-triazol-1-ylmethyl)pentan-3-ol