



Mouse Anti-Clomazone Monoclonal Antibody, clone CMZ (CABT-Z535M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Clomazone-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CMZ
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction Clomazone is an agricultural herbicide, and has been the active ingredient of products named "Command" and "Commence". The molecule consists of a 2-chlorobenzyl group bound to a N-O heterocycle called isoxazolidinone. It is a white solid.

Keywords Clomazone;2-[(2-chlorophenyl)methyl]-4,4-dimethyl-3-isoxazolidinone;Dimethazone
