



# Mouse Anti-Azoxystrobin Monoclonal Antibody, clone AXS (CABT-Z537M)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Azoxystrobin-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	AXS
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

## BACKGROUND

**Introduction**

Azoxystrobin is the ISO common name for an organic compound that is used as a fungicide. It is a broad spectrum systemic active ingredient widely used in agriculture to protect crops from fungal diseases.

---

**Keywords**

Azoxystrobin

---