



Mouse Anti-2,4-Toluene diisocyanate Monoclonal Antibody, clone 2,4-TDI (CABT- Z551M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	2,4-Toluene diisocyanate-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	2,4-TDI
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Toluene diisocyanate (TDI) is an organic compound with the formula $\text{CH}_3\text{C}_6\text{H}_3(\text{NCO})_2$. Two of the six possible isomers are commercially important: 2,4-TDI and 2,6-TDI.

Keywords

2,4-Toluene diisocyanate;Toluene diisocyanate;TDI