



Rabbit Anti-Arabidopsis thaliana Complex I subunit 17.2 kD Polyclonal Antibody (CABT-Z561R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH-conjugated synthetic peptide of 17.2kD Subunit derived from Arabidopsis thaliana AT3G03100.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Unpurified
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:1000-1:2000
Molecular Weight	Expected M.W.: 18 kDa; Apparent M.W.: 17 kDa
Reconstitution	Reconstitution with 150 µl of sterile water.
Cellular Localization	Mitochondria
Format	Lyophilized
Size	150 µg
Preservative	None

Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.
----------------	---

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Complex I is the largest protein complex of the oxidative phosphorylation system in mitochondrial and it catalyzes NADH-quinone oxidoreduction. Complex I represents the main entrance site for electrons into the respiratory electron transfer chain. In Arabidopsis, Complex I have at least 49 subunits and 17.2 kD subunit (AT3G03100) may be one of the subunit.
---------------------	--

Keywords	17.2kD Subunit;B17.2
-----------------	----------------------

GENE INFORMATION

Gene Name	At3g03100
------------------	-----------

Entrez Gene ID	821104
-----------------------	------------------------

UniProt ID	Q01526
-------------------	------------------------