



Mouse Anti-Ganglioside GD3 monoclonal antibody, clone R-24 (CABT-Z60503)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No cross-reaction with other carbohydrate epitopes.
Target	Ganglioside GD3
Immunogen	Purified ganglioside GD3 was used to immunize mice. R-24 hybridoma cells were produced by spleen cell fusion
Isotype	IgG/IgM
Source/Host	Mouse
Species Reactivity	N/A
Clone	R-24
Purification	Serum, Unpurified
Conjugate	Unconjugated
Applications	ELISA,TLC
Molecular Weight	N/A
Format	Liquid
Size	50 µl
Buffer	Dilution: Phosphate buffered saline (pH 7.4) is recommended.
Preservative	None

Storage	-20°C
Ship	Dry ice