



Rabbit Anti-Human RAB12 (Phospho S106) Monoclonal Antibody, clone NKG-S36-10 (CABT-Z618R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Clone	NKG-S36-10
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, Dot Recommended dilution: WB: 1:1000 Dot: 1:1000
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA, pH 7.2.
Preservative	0.01% Sodium azide

Storage	Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Stable for 12 months at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	RAB12 belongs to the small GTPase superfamily and Rab family. RAB12 locates in Golgi apparatus membrane.
Keywords	RAB12;RAB12, member RAS oncogene family;ras-related protein Rab-12;putative Ras-related protein Rab-12

GENE INFORMATION

Gene Name	RAB12
Entrez Gene ID	201475
UniProt ID	Q6IQ22