



Rabbit Anti-Pirimiphos Polyclonal Antibody (CABT-Z647R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	BSA-Pirimiphos-methyl Conjugate
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Caprylic acid extraction purified
Conjugate	Unconjugated
Applications	ELISA Recommended dilution: ELISA: 1/100,000
Format	Liquid
Concentration	Lot specific
Size	0.5 ml
Buffer	Each vial contains 0.5 ml, 0.8mg/ml purified IgG containing 0.02% sodium azide.
Preservative	0.02% sodium azide
Storage	Store at 4 °C or -20 °C. The working antibody solution is stable for at least 7 days at 4 °C.
Ship	Wet ice

BACKGROUND

Introduction	Pirimiphos is an insecticide belonging to the group of organophosphothioates. It is applied in the field on fruits, vegetables, in the storage of grain products and as an additive in paints. It has adverse effects on algae and mammals due to its biological activity as acetylcholine esterase inhibitor.
Keywords	Pirimiphos-methyl;Pirimiphos