



Rabbit Anti-Human PP2D1 Polyclonal Antibody (CABT-Z662R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein Epitope Signature Tag (PrEST) antigen sequence: PSYQMTTDEQQIINSFYTVFREEYAAIEDLFSAINKTEAVRCEYEDTHKAFKAFWRMDR LLGLGRKEVSRVQWSGC
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IHC Recommended dilution: IHC: 1:2500 - 1:5000
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Preservative	0.02% Sodium azide
Storage	Store at 4°C for short term; -20°C for long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Keywords PP2D1;protein phosphatase 2C-like domain containing 1;C3orf48;protein phosphatase 2C-like domain-containing protein 1

GENE INFORMATION

Gene Name	PP2D1
Entrez Gene ID	151649
UniProt ID	A8MPX8
