



Human Anti-Human CD41 (Abciximab) Monoclonal Antibody, clone 7E3 [Biosimilar] (CABT-Z664H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This is a reformatted human IgG1 antibody based on the therapeutic Fab fragment.
Immunogen	Human CD41
Isotype	IgG1
Source/Host	Human
Species Reactivity	Human
Clone	7E3
Purification	Protein A purified
Conjugate	Unconjugated
Applications	BL, FC
Format	Liquid
Concentration	Lot specific
Size	200 µg
Buffer	PBS with 0.02% Proclin 300.
Preservative	None
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.

Ship

Wet ice

BACKGROUND

Introduction CD41 is the integrin alpha chain 2b forming a fibrinogen receptor in heterodimeric complex with beta 3.

Keywords ITGA2B;integrin, alpha 2b (platelet glycoprotein IIb of IIb/IIIa complex, antigen CD41); GT;GTA;CD41;GP2B;HPA3;CD41B;GPIIb;BDPLT2

GENE INFORMATION

Gene Name ITGA2B

Entrez Gene ID [3674](#)

UniProt ID [P08514](#)
