



# Rabbit Anti-Arabidopsis thaliana ALB4 Polyclonal Antibody (CABT-Z707R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	KLH-conjugated synthetic peptide of ALB4 derived from Arabidopsis thaliana AT1G24490.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Purification</b>	Antigen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB Recommended dilution: WB: 1:1000-1:2000
<b>Molecular Weight</b>	55 kDa
<b>Reconstitution</b>	Reconstitution with 150 µl of sterile water.
<b>Cellular Localization</b>	Chloroplast
<b>Format</b>	Lyophilized
<b>Size</b>	150 µg
<b>Preservative</b>	None
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.

## BACKGROUND

**Introduction**

ALB4 is a separate 55 kDa protein and belongs to homologue of the Alb3/Oxa1/YidC family. It is required for proper chloroplast biogenesis.

**Keywords**

ALB4;ALBINA 4;ARABIDOPSIS THALIANA ENVELOPE MEMBRANE  
INTEGRASE;ARTEMIS;ALB3L1.

---

## GENE INFORMATION

**Gene Name**

ALB4

**Entrez Gene ID**

[839065](#)

**UniProt ID**

[Q9FYL3](#)

---