



Human Anti-SARS-CoV-2 NP (N-term) Monoclonal Antibody, clone CO18 (CABT- Z728H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SARS-CoV-2 NP
Isotype	IgG
Source/Host	Human
Species Reactivity	SARS-CoV-2
Clone	CO18
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA (Det) We recommend the following antibodies for ELISA (Capture - Detection): CABT-Z727H - CABT-Z728H
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS or Tris-Gly
Preservative	None
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw

cycles.

Ship

Wet ice

BACKGROUND

Introduction

The 2019 novel coronavirus (SARS-CoV-2), also known as the Wuhan coronavirus, is a contagious virus that causes respiratory infection and has shown evidence of human-to-human transmission, first identified by authorities in Wuhan, Hubei, China, as the cause of the ongoing COVID-19 outbreak. Genomic sequencing has shown that it is a positive-sense, single-stranded RNA coronavirus.

Keywords

SARS-CoV-2;coronavirus;SARS-CoV-2 NP;SARS-CoV-2 Nucleocapsid Protein
