



Rabbit Anti-Arabidopsis thaliana Complex I AT5G08530 Polyclonal Antibody (CABT-Z784R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH-conjugated synthetic peptide derived from N-terminal of AT5G08530 protein in Arabidopsis thaliana AT5G08530.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Unpurified
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:1000-1:2000
Molecular Weight	53 kDa
Reconstitution	Reconstitution with 150 µl of sterile water.
Cellular Localization	Mitochondria
Format	Lyophilized
Size	150 µg
Preservative	None
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after

reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.

Ship	Wet ice
------	---------

BACKGROUND

Introduction	Complex I is the largest protein complex of the oxidative phosphorylation system in mitochondrial and it catalyzes NADH-quinone oxidoreduction. Complex I represents the main entrance site for electrons into the respiratory electron transfer chain. In Arabidopsis, Complex I have at least 49 subunits and 51 kD (AT5G08530) is one of the subunit.
--------------	--

Keywords	AT5G08530;51 KDA SUBUNIT OF COMPLEX I;CI51;NDUFV1;Complex I
----------	---

GENE INFORMATION

Gene Name	At5g08530
-----------	-----------

Entrez Gene ID	830752
----------------	------------------------

UniProt ID	Q9FNN5
------------	------------------------
