



# Rabbit Anti-Arabidopsis thaliana Complex V-1 Polyclonal Antibody (CABT-Z786R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	KLH-conjugated synthetic peptide of Complex V-1 derived from Arabidopsis thaliana ATMG01190, AT2G07698.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Purification</b>	Unpurified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB Recommended dilution: WB: 1:1000-1:2000
<b>Molecular Weight</b>	55 kDa (ATMG01190), 86 kDa (AT2G07698)
<b>Reconstitution</b>	Reconstitution with 150 µl of sterile water.
<b>Cellular Localization</b>	Mitochondria
<b>Format</b>	Lyophilized
<b>Size</b>	150 µg
<b>Preservative</b>	None
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after

reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.

---

Ship	Wet ice
------	---------

---

## BACKGROUND

Introduction	Mitochondrial F0F1-ATP synthase is also called Complex V and it synthesis ATP from ADP and Pi using the proton motive force created by respiratory electron transport. ATP1 (ATMG01190/AT2G07698) is a subunit of mitochondrial F0F1-ATP synthase in Arabidopsis.
--------------	---

---

Keywords	ATP1;ATP SYNTHASE SUBUNIT 1;Complex V-1
----------	---

---

## GENE INFORMATION

Gene Name	ATPA
-----------	------

---

UniProt ID	<a href="#">P92549</a>
------------	------------------------

---