



Rabbit Anti-Arabidopsis thaliana Complex I B7 Polyclonal Antibody (CABT-Z809R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH-conjugated synthetic peptide of Complex I B7 Subunit derived from Arabidopsis thaliana AT2G02050.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:1000-1:2000
Molecular Weight	12 kDa
Reconstitution	Reconstitution with 150 µl of sterile water.
Cellular Localization	Mitochondria
Format	Lyophilized
Size	150 µg
Preservative	None
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after

reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.

Ship	Wet ice
------	---------

BACKGROUND

Introduction	Complex I is the largest protein complex of the oxidative phosphorylation system in mitochondrial and it catalyzes NADH-quinone oxidoreduction. Complex I represents the main entrance site for electrons into the respiratory electron transfer chain. In Arabidopsis, Complex I have at least 49 subunits and B18 (AT2G02050) is one of the subunit.
--------------	--

Keywords	NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 7;At2g02050;Complex I B7
----------	--

GENE INFORMATION

Gene Name	At2g02050
-----------	-----------

Entrez Gene ID	814736
----------------	------------------------

UniProt ID	Q9SKC9
------------	------------------------
