



Goat Anti-PER.C6® HCP Polyclonal Antibody (CABT-Z865G)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PER.C6® HCP
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2°C to 8°C
Ship	Wet ice

BACKGROUND

Introduction	Host cell proteins (HCPs) are process-related protein impurities that are produced by the host organism during biotherapeutic manufacturing and production.
Keywords	HCP;Host Cell Protein