



# Goat Anti-MRC5 HCP Polyclonal Antibody (CABT-Z868G)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	MRC5 HCP
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2°C to 8°C
Ship	Wet ice

## BACKGROUND

**Introduction** Host cell proteins (HCPs) are process-related protein impurities that are produced by the host organism during biotherapeutic manufacturing and production.

---

**Keywords** HCP;Host Cell Protein

---