



Goat Anti-*S. cerevisiae* HCP Polyclonal Antibody (CABT-Z877G)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	<i>S. cerevisiae</i> HCP
Isotype	IgG
Source/Host	Goat
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	Lot specific
Size	500 µl
Buffer	PBS
Preservative	None
Storage	Short Term: 2°C to 8°C
Ship	Wet ice

BACKGROUND

Introduction	Host cell proteins (HCPs) are process-related protein impurities that are produced by the host organism during biotherapeutic manufacturing and production.
---------------------	---

Keywords	HCP;Host Cell Protein
-----------------	-----------------------
