



Rabbit Anti-Human PRADC1 Polyclonal Antibody (CABT-Z882R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein Epitope Signature Tag (PrEST) antigen sequence: DYLYFQVLSPGDIRYIFTATPAKDFGGIFHTRYEQIHLVPAEPPEACGELSNGFFIQDQI ALVERGGCSFLSKTRVVQ
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC Recommended dilution: IHC: 1:20 - 1:50
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	PRADC1, protease-associated domain-containing protein 1, is a novel metabolic-responsive secretory protein that modulates physical activity and adiposity.
Keywords	PRADC1;protease-associated domain containing 1;PAP21;C2orf7;protease-associated domain-containing protein 1;hPAP21

GENE INFORMATION

Gene Name	PRADC1
Entrez Gene ID	84279
UniProt ID	Q9BSG0