



Mouse Anti-Baboon IgG1 Monoclonal Antibody, clone 7H4.G8.G10 (CABT-Z883M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Baboon IgG1
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Baboon
Clone	7H4.G8.G10
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Lot specific
Size	200 µg
Buffer	PBS with 0.02% Proclin 300.
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Immunoglobulin G (IgG) is a type of antibody. IgG is the most common type of antibody found in blood circulation. IgG molecules are created and released by plasma B cells. Each IgG antibody has two paratopes.

Keywords

Ig gamma 1 chain C region;IGHG1;Immunoglobulin heavy constant gamma 1 (G1m marker);IgG1;Immunoglobulin G1;IgG1 hinge region
