



Mouse Anti-Rhesus IgG4 Monoclonal Antibody, clone 2C6 (CABT-Z887M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Rhesus IgG4
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Rhesus
Clone	2C6
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Lot specific
Size	200 μ g
Buffer	PBS with 0.02% Proclin 300.
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Immunoglobulin G (IgG) is a type of antibody. IgG is the most common type of antibody found in blood circulation. IgG molecules are created and released by plasma B cells. Each IgG antibody has two paratopes.

Keywords

IGHG4;immunoglobulin heavy constant gamma 4 (G4m marker);IGG4;immunoglobulin heavy constant gamma 4
